

NAME

dicom2scn – medical image conversion

SYNOPSIS

```
dicom2scn [ -b frame ] [ -e frame ] [ -r ]  
          [ -i ] inputfile [ outputfile ]
```

```
dicom2scn -h
```

```
dicom2scn --help
```

DESCRIPTION

dicom2scn converts sequences of 2D DICOM image files to the SCN volume format (see **scn(5)**). It was developed to work with volumetric T1 DICOM exports from the Elscint MR scanner at Unicamp's Hospital. It is quite improbable that it will convert DICOM files from elsewhere.

DICOM NAMING

All DICOM slices belonging to a study should be in a same directory, named with a fixed prefix and a numeric suffix, e.g.:

```
SOME_PATIENT1  
SOME_PATIENT2  
(...)  
SOME_PATIENT162
```

The *inputfile* parameter can be any slice filename (e.g.: */some/directory/SOME_PATIENT85*) or the path to a directory that contains only DICOM files. If no *outputfile* is given, the conversion output (if any) is written to *dicom.scn*

OPTIONS

```
-b frame  
    begins the conversion attempt at the frame-th slice.  
  
-e frame  
    ends the conversion attempt at the frame-th slice.  
  
-r  
    writes a conversion report to dicom.txt  
  
-i  
    lists the image sequences on these DICOM files and some of their properties.  
  
-h, --help  
    prints a brief usage help to the console and exits.
```

LIMITATIONS

It may be common to have multi-study exports, with different studies in the same file sequence. With no **-b** or **-e** parameters, dicom2scn will fail with a quite verbose error message if it crosses the boundary between two studies. The error message will provide suggestions for **-b** and/or **-e** study boundaries. dicom2scn requires that slices be no wider than 2048 pixels.

AUTHORS

dicom2scn was written by Felipe Bergo <bergo@seul.org>. Information about IVS and associated utilities can be found at

<http://www.ic.unicamp.br/~afalcao/ivs>

SEE ALSO

ivs(1), ana2scn(1), scn2ana(1), scntool(1), scn(5)